

Collars

A. Nylon or Leather Collars

To determine the proper size, measure the circumference of your pet's neck and add 2". Collars should fit high on pet's neck without fitting too loosely.


$\frac{1}{2}$ "- $\frac{3}{4}$ " width is best for puppies and small dogs. $\frac{3}{4}$ " width is best for small to medium dogs. 1"-1 $\frac{1}{2}$ " width is best for large dogs.

C. Choke Chain Training Collars

To determine the proper size, measure the circumference of your pet's neck and add 2"-3". The effectiveness of the correction will be diminished if collar is too loose.


B. Prong Chain Training Collars

To determine the proper size, measure around the dog's neck and add 1". Prong Collars should fit snug and high on your dog's neck. Millimeter (mm) size refers to the thickness of the metal. For a longer collar in the thickness that you prefer, add extra links.



Beds

When dog is laying down, measure pet's body length from base of tail to tip of nose, and add a minimum of 9-12". For bolster or cuddler style beds, add a minimum of 5-7".



Selecting a crate mat

To determine the proper size crate mat, refer to the base dimensions of your crate or kennel. If dimensions do not match exactly, we recommend selecting a crate mat $\frac{1}{2}$ "-1" smaller than base dimensions of your crate.

Crates

To determine proper crate size, measure dog's height to the shoulders, and length from chest to base of tail, then add 4"-6" to each measurement. See below for estimations by breed.

XS (x-small)

Cat, Yorkshire Terrier, Chihuahua, Maltese, Italian Greyhound

S (small)

Beagle, Boston Terrier, Cocker Spaniel, Pug

M (medium)

Basset Hound, Bulldog, Keeshond, Norwegian Elkhound, Dalmatian

ML (med-lg)

Boxer, Standard Poodle, Golden Retriever, Labrador


L (large)

Afghan, Doberman, German Shepherd, Irish Setter, Rottweiler, Collie

XL (x-large)

Great Dane, Great Pyrenees, Greyhound, Mastiff, Saint Bernard

The sizes indicated above are recommended for travel. For household use, select a size or two larger. This will allow the dog ample room to stand, turn around and lie down. If your dog's breed is not listed above, compare the dog's body size to a similar-sized breed and follow measuring instructions.




Apparel

A. Coats, Sweaters and Shirts

Measure length of pet's back from the base of the neck to the base of the tail. If the measurement is between sizes, we recommend going up one size. Also, if your dog has a stout build or is overweight, we recommend going up a minimum of one size for proper fit.

B. Harnesses and Vests

Measure girth around the broadest part of the chest and add 2".



Muzzles

A. Basket

Measure midway between tip of nose and eyes. Add 2" for small/medium breeds and 3" for large breeds.

B. Open ended fabric muzzles

Measure just above the nose.

Pet Doors

Use the chart below to select the proper size pet door based on the dog's adult body weight. If the dog is at or slightly heavier than the weight indicated, select the next larger size. For puppies, use the breed list to select according to the full-grown size of the breed. If the dog is overweight, measure across the widest part of the chest and select a door with a flap opening that is slightly wider than the width measurements indicated below.

		Pet Safe Door Flap Opening	
		W x H	
Small	For pets up to 10 lbs	4 $\frac{1}{2}$ "W	7"H
	Cat, Chihuahua, Toy Poodle		
Medium	35 lbs	8"W	11 $\frac{1}{2}$ "H
	Terrier, Spaniel, Beagle and similar-sized breeds		
Large	75 lbs	10 $\frac{1}{2}$ "W	15 $\frac{1}{2}$ "H
	Shepherd, Boxer, Collie and similar-sized breeds		
Extra Large	120 lbs	13"W	18 $\frac{1}{2}$ "H
	Great Dane, St. Bernard and similar-sized breeds		